

Intestinal organoid cultures

 Zuri Sullivan  Ruslan Medzhitov

Updated date: May 18, 2021

 An abbreviated version of this protocol was published in Science in Mar 2021

$\gamma\delta$ T cells regulate the intestinal response to nutrient sensing

DOI: 10.1126/science.aba8310

Detailed protocol

Adapted from Sato et al *Nature* 2009 and Sato et al *Gastroenterology* 2011

Basal organoid medium

500mL advanced DMEM/F12
100U/mL Penicillin/Streptomycin (100X)
10mM HEPES (100X)
200mM GlutaMAX (100X)
N2 supplement (100X)
B27 supplement (50X)
1mM N-acetylcysteine (500X 81.6mg/mL stock in dH₂O)

Store at 4 degrees for ~4 weeks

ENR medium

50mL basal medium
50ng/mL rmEGF (10,000X 500ug/mL stock in 0.1% BSA/PBS)
100ng/mL rmNoggin (1000X 100ug/mL stock in 0.1% BSA/PBS)
1ug/mL rhRspo1 (1000X 1mg/mL stock in 0.1% BSA/PBS)

Store at 4 degrees for 1 week

Crypt isolation

1. Isolate 20cm of small intestine and open longitudinally, rinse in cold PBS
2. Chop tissue into 2-5mm pieces and transfer to 25mL 2mM EDTA/PBS
 - a. Incubate 30 minutes on ice
 - b. During incubation, prepare ENR medium
 - c. Prewarm 24-well non-TC treated plate
3. Remove EDTA medium with 10mL pipette
4. Resuspend tissues vigorously in 20mL cold 0.1%BSA/PBS using 10mL pipette for >30s
 - a. Let settle and remove supernatant (villous fraction)
5. Resuspend pellet vigorously in 10mL cold 0.1%BSA/PBS using 10mL pipette and shake for 30s. Pass through 70uM strainer
6. Centrifuge at 200x g for 3 minutes at 4 degrees and gently discard supernatant
 - a. Prewarm 24-well non-TC treated plate at 37 degrees
 - b. Place Matrigel on ice
7. Resuspend crypts in 10mL cold basal medium, transfer to 15mL conical, and place on ice
 - a. count crypts using hemocytometer
 - b. aim to plate ~500 crypts per well
8. Centrifuge 200x g for 3 minutes at 4 degrees and remove supernatant
9. Resuspend crypts at a concentration of 500 crypts/100uL in 50% ENR medium/50% Matrigel
10. Quickly plate 50uL/well in pre-warmed 24-well non-TC treated plate
11. Incubate at 37 degrees for 10 minutes to set Matrigel
12. Add 500uL ENR medium and change medium 3x per week

Organoid maintenance

Passage organoids 1:5 once per week

1. Remove culture medium and add 500uL basal medium
 - a. Thaw Matrigel on ice and pre-warm 24-well non-TC treated plate
2. Use 1000uL pipette to mechanically disrupt organoids + matrigel and transfer to 15mL falcon tube
 - a. Further dissociate organoids using fire-polished Pasteur pipette
3. Add 10mL basal medium and centrifuge at 200x g for 2 minutes
4. Discard supernatant and resuspend in 125uL ENR medium + 125uL Matrigel
5. Quickly plate 50uL/well in pre-warmed 24-well non-TC treated plate
6. Incubate at 37 degrees for 10 minutes to set Matrigel
7. Add 500uL ENR medium and change medium 3x per week

How to cite:(Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Sullivan, Z. and Medzhitov, R. (2021). Intestinal organoid cultures. Bio-protocol Preprint. [bio-protocol.org/prep1098](https://doi.org/10.21203/rs.3.rs-1098).

2. Sullivan, Z. A., Khoury-Hanold, W., Lim, J., Smillie, C., Biton, M., Reis, B. S., Zwick, R. K., Pope, S. D., Israni-Winger, K., Parsa, R., Philip, N. H., Rashed, S., Palm, N., Wang, A., Mucida, D., Regev, A. and Medzhitov, R.(2021). $\gamma\delta$ T cells regulate the intestinal response to nutrient sensing. Science 371(6535). DOI: [10.1126/science.aba8310](https://doi.org/10.1126/science.aba8310)

Copyright: Content may be subjected to copyright.